

# COVID-19 X-RAY DETECTION

Vijay Saddi  
August 2020



# AGENDA



USE CASE



DATA SET



DATA PREP



MODEL BUILDING



EVALUATION






# COVID-19 BACKGROUND

- The 2019 outbreak of Covid-19 has kept the world on lockdown with the effort to stop the spread
- It is a highly contagious illness which affects the lungs of the respiratory system
- Transmitted via close proximity via small droplets produced by coughing , sneezing & talking
- X-rays and CT scans of Covid-19 victims can depict the damage inflicted on the lungs



# RADIOLOGY



- Chest x-rays are useful in the follow-up of the disease (it is insensitive early on)
- The image on the left is an 83-year old male with Covid-19
- Image findings:
  - Admission: Ill defined bilateral alveolar consolidation
  - 4-hrs: Radiology worsening, affected lower lobes
  - 24-hrs: Bilateral alveolar consolidation<sup>1</sup>
  - 48-hrs: Further radiology worsening
  - 72-hrs: Panlobar affectation<sup>2</sup>, with typical radiology findings of ARDS<sup>3</sup>; 24 hrs later patient passed away

1. Infers an alveolar spread of disease, commonly due to pneumonia  
2. Alveolar destruction within the lobule  
3. Acute Respiratory Distress Syndrome: fluid buildup in the alveoli



# CURRENT STATISTICS ON COVID-19

Last updated: August 08, 2020, 20:24 GMT

Coronavirus Cases:

**19,725,620**

Deaths:

**727,540**

Recovered:

**12,659,896**

- High rate of infection
- High rate of recovery (95%)
- 5% death rate (August 2020)

## ACTIVE CASES

**6,338,184**

Currently Infected Patients

**6,273,265** (99%)  
in Mild Condition

**64,919** (1%)  
Serious or Critical

[Show Graph](#)

## CLOSED CASES

**13,387,436**

Cases which had an outcome:

**12,659,896** (95%)  
Recovered / Discharged

**727,540** (5%)  
Deaths

[Show Graph](#)




# QUESTION

With the emergence of Deep Learning for AI in healthcare, can we build a model to help identify Covid-19 by evaluating lung x-rays?



# DATA SET




- 2 main data sources for Covid-19 and Normal x-ray images
- Covid-19 x-rays were taken from Github:  
<https://github.com/ieee8023/covid-chestxray-dataset>
- Normal x-ray images were taken from Kaggle:  
<https://www.kaggle.com/paultimothymooney/chest-xray-pneumonia>
- There were 180 images per group, with a total of **360**





# EXPLORATION & QUALITY ASSESSMENT

- Image file sizes range from 400-600kb each
- Images are of high quality (720p-1080p resolution)
- Each image has different dimensions (width vs height)
- RGB & greyscale formats






# DATA PRE-PROCESSING

- Github images were filtered to select only Covid-19 diagnosed, and Posteroanterior view
- From the entire pool of images, the *80/20 split* was implemented to divide the training/test data
- **Training:** 144 images per category (Covid & Normal)
- **Testing:** 36 images per category (Covid & Normal)

```
df_metadata[(df_metadata["finding"]=="COVID-19") & (df_metadata["view"]=="PA")]
```

	patientid	offset	sex	age	finding	RT_PCR_positive	survival	intubated	intubation_prese
0	2	0.0	M	65.0	COVID-19	Y	Y	N	N



# FEATURE ENGINEERING

- Prior to model training, augmentation methods were required to account for variations in the image makeup:
  - **Normalization:** images rescaled by 1/255 to ensure pixel values range between 0 and 1
  - **Shear range:** shear angle for variations in image rotation
  - **Zoom range:** account for variations in zoom between images
  - **Horizontal flip:** reverses the rows/columns of pixels in random images – as lung images show symmetry in nature
  - **Target size:** resize all images to standard 224x224

```
# prepping & augmenting the images for loading into the model as inputs for training
train_data = image.ImageDataGenerator(
    rescale = 1./255, #normalize images
    shear_range = 0.2,
    zoom_range = 0.2,
    horizontal_flip = True
)
```

```
# converting image size, set the batch size (standard 224 by 224)
# feed the augmented images from previous step into the generator
train_generator = train_data.flow_from_directory(
    'Covid_Dataset/Train',
    target_size = (224,224),
    batch_size = BS,
    class_mode = 'binary'
)
```



# DEEP LEARNING MODEL

- A Convolutional Neural Network (CNN) is best suited for image processing
- Summary:
  - 4 convolutional layers
  - 3 pooling layers (2,2)
  - 4 dropout layers
- Dropout Layer: reduction by 25% to prevent overfitting in first 3 layers
- Loss function: binary cross-entropy to classify between 2 classes
- Optimizer: Adam shows better performance on deep learning models

```
model = Sequential()

model.add(Conv2D(32, kernel_size=(3,3), activation='relu', input_shape=(224,224,3)))
model.add(Conv2D(64, (3,3), activation='relu'))
model.add(MaxPool2D(pool_size=(2,2)))
model.add(Dropout(0.25))

model.add(Conv2D(64, (3,3), activation='relu'))
model.add(MaxPool2D(pool_size=(2,2)))
model.add(Dropout(0.25))

model.add(Conv2D(128, (3,3), activation='relu'))
model.add(MaxPool2D(pool_size=(2,2)))
model.add(Dropout(0.25))

model.add(Flatten())
model.add(Dense(64, activation='relu'))
model.add(Dropout(0.5))
model.add(Dense(1, activation='sigmoid'))

model.compile(loss=keras.losses.binary_crossentropy, optimizer='adam', metrics=['accuracy'])
```




Model: "sequential\_3"

Layer (type)	Output Shape	Param #
conv2d_9 (Conv2D)	(None, 222, 222, 32)	896
conv2d_10 (Conv2D)	(None, 220, 220, 64)	18496
max_pooling2d_7 (MaxPooling2)	(None, 110, 110, 64)	0
dropout_9 (Dropout)	(None, 110, 110, 64)	0
conv2d_11 (Conv2D)	(None, 108, 108, 64)	36928
max_pooling2d_8 (MaxPooling2)	(None, 54, 54, 64)	0
dropout_10 (Dropout)	(None, 54, 54, 64)	0
conv2d_12 (Conv2D)	(None, 52, 52, 128)	73856
max_pooling2d_9 (MaxPooling2)	(None, 26, 26, 128)	0
dropout_11 (Dropout)	(None, 26, 26, 128)	0
flatten_3 (Flatten)	(None, 86528)	0
dense_5 (Dense)	(None, 64)	5537856
dropout_12 (Dropout)	(None, 64)	0
dense_6 (Dense)	(None, 1)	65
Total params: 5,668,097		
Trainable params: 5,668,097		
Non-trainable params: 0		



# MODEL PERFORMANCE

- Based on the confusion matrix, the model has a high accuracy when predicting Covid-19 and Normal x-rays
- **97%** accuracy when predicting Covid-19 x-rays
- **94%** accuracy when predicting normal x-rays



# PERFORMANCE INDICATORS

- **F1-Score** is used as the best performance indicator for binary classification problems
- **Precision** – correctly identified 97% of positive cases from all the predicted positive cases
- **Recall** – correctly identified 94% of positive cases from all the actual positive cases

	precision	recall	f1-score	support
Covid	0.97	0.95	0.96	37
Normal	0.94	0.97	0.96	35
accuracy			0.96	72
macro avg	0.96	0.96	0.96	72
weighted avg	0.96	0.96	0.96	72
Precision score: 0.9714				
Recall score: 0.9444				
<u>F1 Score: 0.9577</u>				



# MODEL EVALUATION




96%

The CNN model yields an F1-score  
of 96% on test data





# TRAINING HISTORY



- Does the model overfit to the training data?
- Based on the plot history, the accuracy and loss follows the same trend, although does not overfit!
- Throughout each Epoch, there are slight variations in loss and accuracy

# SUMMARY/NEXT STEPS



## CNN MODEL

- Successful working model
  - 96% Accuracy
  - Does not overfit



## LIMITATION

Data set is very small for large scale use; more data needed



## NEXT STEP

- Using ImageNet to import a larger data set for training
- Compare model with a transfer learning model results



THANKS!



# CONTACT

Does anyone have any questions?

Email: [vijay.saddi@outlook.com](mailto:vijay.saddi@outlook.com)

LinkedIn: <https://www.linkedin.com/in/vijay-ss/>

Github: <https://github.com/vijay-ss>

Portfolio: [https://vijay-ss.github.io/Vijay\\_Portfolio/](https://vijay-ss.github.io/Vijay_Portfolio/)

